

Document No: GENE-WM391	Version No: 2
	Issue date: 30/09/2022



OSTEOGENESIS IMPERFECTA GENE PANEL - TECHNICAL INFORMATION

Design: The Osteogenesis imperfecta (OI) gene panel was designed as part of a custom probe set from Twist Bioscience (TE-98175847) to cover 19 genes associated with OI. This panel design provides coverage of 100% of the target coding regions and flanking intronic sequences (+/- 15bp) for the 18 genes listed below.

Method: Library preparation and target enrichment was performed using the custom designed TE-98175847 probe set (Twist Bioscience) and Nextera Flex for Enrichment (Illumina). Sequencing was performed using a 150bp paired-end sequencing kit on a MiSeq (Illumina). All stages of the workflow were performed according to the manufacturer's instructions.

Coverage criteria: For each sample reported, >95% of the target coding and flanking intronic regions of the 19 genes was covered to a minimum depth of 20 reads (20x). Any regions of the COL1A1 and COL1A2 genes not covered to 20x depth were flagged for follow-up Sanger sequencing. Specific details of coverage and depth for individual tests are available from the laboratory on request.

Variant identification and interpretation: Sequence data were mapped and variants identified using GenomeAnalysisToolKit (GATK) and NextGENe (Softgenetics) with hg19 (GRCh37) human genome as the reference. Variants identified were subsequently classified according to recent ACGS Best Practice Guidelines for Variant Classification using all available evidence. Any clinically significant variants were confirmed by Sanger sequencing.

Variant reporting: Variants were reported according to HGVS guidelines using the accession numbers listed below. Variants categorised as non-pathogenic or likely non-pathogenic were filtered out at the variant interpretation stage and were not included in the clinical report. Details of these variants are available from the laboratory on request.

Genes included and associated sequence accession numbers:

BMP1 (NM_006129.4); COL1A1 (NM_000088.3); COL1A2 (NM_000089.3); CREB3L1 (NM_052854.3); CRTAP (NM_006371.4); FAM46A (TENT5A) (NM_017633.2); FKBP10 (NM_021939.3); IFITM5 (NM_001025295.2); KDELR2 (NM_006854.4); P3H1 (NM_022356.3); PLOD2 (NM_182943.2); PLS3 (NM_005032.6); PPIB (NM_000942.4); SERPINF1 (NM_002615.6); SERPINH1 (NM_001235.3); SP7 (NM_152860.1); SPARC (NM_003118.3); TMEM38B (NM_018112.2); WNT1 (NM_005430.3).

Authority For Issue: Christine Black	Page 1 of 1
Document printed from Q-Pulse on 04/10/2022 by Black, Holly	

This is a controlled document: This copy is valid on day of print only, after which the user must ensure that this is the correct version by comparing against the current document details in Q-Pulse.