

How to use Bluetooth & app?

Pair hearing aids to your smartphone

1. Go to your phone's home screen and tap on the “**Settings**” icon.
2. In the settings menu, tap on the **Bluetooth** icon.
3. Turn Bluetooth **On**.
4. Switch the hearing aids off and on.
5. After several seconds, the name of your Phonak hearing aids (L-NHS hearing aid / R-NHS hearing aid) will appear on your phone screen. Tap on it to select.
 - a. If you have two hearing aids, only the name of one hearing aid will appear. This is normal. Tap on it to select and it will pair to both hearing aids.
6. You will now see a message, “Bluetooth Pairing Request” Click “Pair”.
7. A connection **beep** will be heard in the paired hearing aid.

Download myPhonak app

You can download the **myPhonak** app from Apple iTunes or Google Android app stores.

Type in “myPhonak”. The app you are looking for will be green in color. This is a **free app**. Walk through the welcome screens and you will have to read and accept the privacy notice. Ensure that your Bluetooth is ON.

You may also download the app from the QR code below:



Video Guide

Pair hearing aids to Android phone.



Pair hearing aids to iPhone.



More support

For more information, visit Phonak's support page:
<https://bluetooth.phonak.com>

